

# ZShape

*New generation*  
of body shaping



**Zimmer**  
Aesthetic Division

*Our new generation*

## **ZSHAPE: DELICATE WRINKLE CARE AND BODY CONTOURING**

ZShape is a next-generation radiofrequency device designed to deliver superior results with unparalleled patient comfort.

### **EXPERIENCE SKIN TIGHTENING WITH DEEP HEAT**

RF energy effectively targets deep skin layers, prompting collagen fibers to contract, visibly rejuvenating and revitalizing the skin.





## VACUUM

Enhanced blood flow and lymph drainage

The medium and large handpieces incorporate a vacuum function that stimulates metabolism, enhances blood circulation, and promotes lymphatic drainage effects.

## RADIO FREQUENCY (RF)

Enhanced skin tone and reduced cellulite

Surface and vacuum chamber electrodes ensure uniform delivery of RF energy across the entire treatment area, and to the optimal depth.

The synergy of RF and Vacuum technologies boosts treatment success and patient satisfaction.



# HANDPIECES

PERFORM AN RF TREATMENT LIKE NO OTHER

## MULTIFUNCTIONAL HANDPIECES

ZShape's unique handpieces have each been specifically designed to accommodate a range of treatment areas.

*small*

*medium*

*large*



## INTELLIGENT FEATURES FOR EFFORTLESS TREATMENT SESSIONS

### THE HIGHEST LEVEL OF FUNCTIONALITY FOR ENHANCED COMFORT AND SAFETY

Each handpiece is equipped with an LED ring light that provides real-time feedback on the therapy.

#### COUPLING INDICATOR:

If optimal skin contact is lost the handpiece begins to flash indicating that the position needs adjustment or contact gel is required.

#### TEMPERATURE INDICATOR:

Real-time temperature feedback ensures safety by preventing over or under heating of the treatment area.

blue LED = temperature is too low

green LED = optimal treatment temperature

orange LED = increased temperature

In addition, the precise temperature is continuously displayed on-screen.

#### CONTROL AT YOUR FINGERTIPS:

Seamlessly control the therapy by adjusting RF energy intensity directly from the handpiece, eliminating the need to adjust settings on the device.





# PATIENT AND PRACTITIONER BENEFITS



## IMMEDIATE RESULTS

The positive effects of reducing fine lines and wrinkles are typically visible to patients after their first treatment.



## PAIN-FREE

ZShape is a very comfortable treatment, even for people sensitive to pain. Patients feel a pleasant warm sensation and a gentle massaging effect from the vacuum.



## NON-INVASIVE

Unlike many other treatments, ZShape is a completely non-invasive treatment method. No needles, just electrodes that glide gently over the skin.



## NO DOWNTIME

Patients can immediately return to work or attend their evening event directly after the treatment.



## MINIMAL TIME COMMITMENT

Prior to therapy, a brief skin cleansing is all that's needed. No post-treatment care is required, reducing treatment time for both staff and patients.



## USER INTERFACE

- » Intuitive design
- » Presets for different body regions
- » Individual adjustment options for maximum efficiency
- » Space to save favorites

Simple operation via touch or with the rotary knob. Select the desired treatment region with a tap of your finger. The appropriate handpiece and the settings needed to start the therapy are then pre-selected.

In addition to the presets, individual adjustments can be made at any time and saved under 'Favorites'.

## PROVIDE PERSONALISED TREATMENTS

- » With 1 MHz, 2 MHz and 4 MHz, ZShape offers 3 different frequencies for effective treatment.
- » Constant vacuum or choose a range of one of 3 different pulse modes.
- » 10 vacuum levels.



# APPLICATION AREAS

## *ZShape*

- TARGETED SMOOTHING OF WRINKLES AND FINE LINES
- BODY CONTOURING IN AREAS SUCH AS THE ABDOMEN, BUTTOCKS, ARMS, THIGHS AND CALVES.
- SKIN TIGHTENING E.G. AFTER PREGNANCY AND WEIGHT LOSS
- TIGHTENING OF SAGGY ARMS
- REDUCTION OF CELLULITE



DÉCOLLETÉ  
NECK



ARMS

ABDOMEN

BUTTOCKS

THIGHS

CALVES

**Small handpiece**

**Medium handpiece**

**Large handpiece**



# COMBINATIONS WITH ZSHAPE

## COMBINE THE DIFFERENT MODES OF OUR BODY SOLUTION SYSTEMS FOR BEST RESULTS!

### COMBINATION 1

If the goal is both muscle growth and body contouring, Zimmer's cutting-edge ZShape and ZStim devices make a perfect pair.

A combination of both treatments in one session, two times a week for a total of 10-12 treatments over a period of 6-8 weeks is recommended.



### COMBINATION 2

ZLipo, ZWave, ZShape and ZStim, a perfectly coordinated team for your practice. Fat reduction with ZLipo, rapid removal of fat tissue with ZWave, body shaping with skin-tightening and muscle-building effects with ZShape and ZStim.

We recommend carrying out the first treatment approx. 1 week after the last ZLipo treatment together with ZWave. From the second week onwards, we recommend treatments with ZShape and ZStim. Every week for at least 8-10 weeks.





## ZIMMER AESTHETIC DIVISION - A COMPLETE BODY SOLUTION

HIGHLY INNOVATIVE  
CRYOLIPOLYSIS™  
THERAPY



ZLIPO



ZWAVE

RADIAL  
SHOCKWAVE  
THERAPY



ZSTIM

ELECTROTHERAPY



# ZShape

*New generation*  
of body shaping

## TECHNICAL DATA

**LCD touchscreen**

**3** handpieces

**4** different pulse modes

**0–10** vacuum level

**1 MHz, 2 MHz, 4 MHz** output frequency

**1–20 RF intensity** levels

**118<sub>H</sub> 57<sub>W</sub> 56<sub>L</sub>** Dimensions (cm)

**30<sub>kg</sub>** Weight

**Zimmer**  
AestheticDivision

follow us



Visit our homepage  
[www.zimmer-aesthetics.de](http://www.zimmer-aesthetics.de)



Zimmer MedizinSysteme GmbH  
Junkersstraße 9  
89231 Neu-Ulm, Germany  
Tel. +49 7 31. 97 61-0  
Fax +49 7 31. 97 61-118  
[www.zimmer-aesthetics.de](http://www.zimmer-aesthetics.de)

EN 10105545 non medical | Version 2 | 1124 | Subject to change |  
Information on possible risks, undesirable side effects, and warnings can be found in the user manual.

